

IK AND

rers of Biochemistry Analyzers,

KANAD VIDYU

Manufacturers of Biochemistry Analyzers and Colorimeters



AMEDI A

DIGITAL HAEMOGLOBINOMETERS





DRY BATH INCUBATOR



HEMA

- Haemoglobinometer with single memory for Hb factor, with CALIBRATION facility.
- Solid State (LEDs) Life Long light source

HEMA 2

- Haemoglobinometer with dual memory for Hb factor as well as Glucose factor with CAL facility
- Solid State (LEDs) Life Long light source

DB 18

- Dry Bath Incubator with 18 Cuvettes/test tubes Capacity
- •Square (12.5mm X 12.5mm) & Round (Dia. 12mm) both types accommodated
- Microprocessor controlled precise Temp up to 37°C.

HANS DIGITAL COLORIMETERS





PRODUCT NAME	WAVELENGTHS	PRECISION & MODES
HANS 251	Five Wavelengths 405, 505, 530, 580, 630	Up to 2 Decimal Precision (0.00), Absorbance Mode Only
HANS 253	Five Wavelengths 405, 505, 530, 580, 630	Up to 2 Decimal Precision (0.00) Absorbance, Concentration & Factor Modes
HANS 291	Nine Wavelengths 405, 430, 465, 505, 530, 565, 585, 630, 700	Up to 2 Decimal Precision (0.00), Absorbance Mode Only
HANS 293	Nine Wavelengths 405, 430, 465, 505, 530, 565, 585, 630, 700	Up to 2 Decimal Precision (0.00) Absorbance, Concentration & Factor Modes
HANS 391	Nine Wavelengths 405, 430, 465, 505, 530, 565, 585, 630, 700	Up to 3 Decimal Precision (0.000), Absorbance Mode Only
HANS 393	Nine Wavelengths 405, 430, 465, 505, 530, 565, 585, 630, 700	Up to 3 Decimal Precision (0.000), Absorbance, Concentration & Factor Modes

KANAD VIDYUU

Manufacturers of Biochemistry Analyzers and Colorimeters



BIOCHEMISTRY mini-ANALYZER

(Programmable Advanced Colorimeter)





SHREERANG 11

- >VIS-5 wavelengths 450, 505, 530, 580, 630
- ➤ Semi-Auto Analyser 35 Test program memory
- ➤ End Point and Fixed Time modes
- >Zero Maintenance- Solid State (LEDs) Life Long Sources
- >Performs Electrolytic (e.g Ca, Mg, Na, K, Cl) and
- > Turbidimetric (RA, ASO, CRP, HbA1c) efficiently.

KANAD PHOTON 35

- >VIS-9 wavelengths 405,430,465,505, 530, 565, 630,700
- ➤ Semi-Auto Analyser 35 Test program memory
- ➤ End Point and Fixed Time modes
- **≻Zero Maintenance- Solid State (LEDs) Life Long Sources**
- > Performs Electrolytic (e.g Ca, Mg, Na, K, Cl) and
- > Turbidimetric (RA, ASO, CRP, HbA1c) efficiently.

SEMI-AUTO BIOCHEMISTRY ANALYZER

(COMPLETE SOLUTION FOR CLINICAL CHEMISTRTY)







MULTI-POINT CALIBRATION,

TURBIDIMETRY



MODEL NAME	KANAD PHOTON 6P	KANAD PHOTON 7P	KANAD PHOTON 8MPC	KANAD PHOTON 9 MPC
MEMORY INPUT	25 Pre-programed tests & 25 user defined programs	100 user defined programs.	150 User defined Programs	200 user defined programs
WAVELGENGTHS (LED source)	340, 405, 465, 505, 545, 630	340, 405, 450 , 505, 545, 580, 630	340, 405, 450, 505, 540, 580, 630, 700	340, 405, 450, 470, 505, 540, 580, 630, 700
ANALYSIS	KINETIC, FIXED TIM	E, END POINT,	KINETIC, FIXED TIME, END POINT,	

COMMON FEATURES

MODES

REAGENT VOLUME
 SIZE
 : OPEN SYSTEM With ONLY 250 μI (0.250ml) reagent volume.
 : 260mm X 240mm X 96mm (L X W X H) (Less than 1ft X 1 ft)

• INPUT / OUTPUT INTERFACE: SPILL PROOF Keyboard, Serial port (RS 232)

LINEARITY : 0.000 - 3.500 Ab
TEMPERATURE : RT and 37°C
DRIFT : 0.001/hr
BATTERY BACKUP : 10 hrs

TURBIDIMETRY



KANAD VIDYUU

Manufacturers of Biochemistry Analyzers and Colorimeters

ADVANTAGES OF CUVETTE MODE ANALYZER OVER ASPIRATION MODE

- No Bubble Formation Errors
- **☞ No Pump calibration**
- No tube Change
- Very low reagent volume reduces expenditure
- No Moving part hence no Mechanical failure
- Very fast readings as No Cleaning Cycle, No Mechanical Errors, No Reagent wastage, No warm-up period
- More number of tests performed including Hb, turbidity, Electrolytes & ELISA readings
- No Vibrations of moving parts
- Noiseless operation

@ KANAD VIDYUT & Technology

- $\ensuremath{\mathfrak{S}}$ Solid state light sources (or LEDs) have a huge life of more than 1,00.000 hrs as compared to halogen lamp which has only 4000hrs life.
- As we use Solid state light source, we don't require filters to get the required wavelengths for measurements.
- © Solid state light sources (LEDs) emits light at same wavelengths and intensity throughout its life. This reduces maintenance
- TAS Halogen lamp is not used, heating problems are eliminated.
- Tuniquely designed optical uses no moving parts in the machine while giving better performance than best.
- Full metal powder coated Aluminum Box gives high sturdiness and long life to instrument.
- Tuniquely designed software programs help reduce errors while performing operations.







Manufactured By :: KANAD VIDYUT

Office: 10, Deenanath Mangeshkar Hospital Road, Erandawane Pune-411004

Contact No:: (+91)(India)-9850848722, 020-25462809

Email: kanad.pune@gmail.com

Website:: www.kanadphoton.com, www.kanadvidyut.com

